

SPLIT II

STEEL

/FRAME	Steel, powder coating
/CONTROLLER	PG/China
/MOTOR	250W×2, brush motor
/MAX SPEED	7 KM/H
/TRAVEL MILEAGE	15 KM
/BATTERY	Lead acid 2×12V, 20Ah
/FRONT WHEEL	8"×2", PU tire
/REAR WHEEL	12.5" pneumatic/PU, steel spoke
/ARMREST	Flip back slope armrest, PU armpad
/FOOTREST	Detachable footrest with heel support
/SEAT	Tension nylon with sponge, thick cushion
/SPECIAL	Electromagnetic brake; Safety belt; Backrest half folding



ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SPLIT II	15"	108	63	97	46	38	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81
	18"	108	63	97	46	45	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81
	20"	108	63	97	46	50	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81