

SPLIT

ALUMINUM

/ FRAME	Aluminum, powder coating
/ CONTROLLER	PG/ China
/ MOTOR	200W×2, brush motor
/ MAX SPEED	7 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lithium 24V 20Ah (Lead acid 20Ah battery optional)
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12.5" pneumatic/PU, aluminum alloy
/ ARMREST	Flip back slope armrest, PU armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Breathable backrest and cushion
/ PACKAGE	95pcs/20"GP,245pcs/40"HQ
/ OPTION	Height adjustable armrest; Color upholstery
/ SPECIAL	Backrest half folding; Safety belt; Electromagnetic brake



ITEM NO.	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SPLIT	108	63	97	34	45	44	44	53	135	8"	12.5"	40,5/33,5	87x38x81