



INVACARE® PERFECTO₂™ V
OXYGEN CONCENTRATOR



Yes, you can.®

The **Invacare® Perfecto2™ V oxygen concentrator** is an economical addition to the premier Invacare stationary concentrator line. This 5-liter concentrator is small, quiet and lightweight and comes standard with the **Invacare® SensO2® Oxygen Sensor**. The Perfecto2 V unit incorporates many of the best features of the Invacare® Perfecto2™ Concentrator, including its contemporary, patient preferred external design and color.

Features

Self diagnostic electronics minimize troubleshooting time.

Top handle for ease of mobility.

Battery-free power loss alarm, compressor 35 psi pressure-relief valve and low flow alarm.

Provider Benefits

Robust design comes standard with Invacare® SensO2® oxygen monitor to simplify preventative maintenance schedules and lower costs.

All Perfecto2 V stationary concentrators include a HomeFill system port for inventory flexibility.

Patient-accessible reset button helps minimize emergency service calls.

Patient Benefits

User friendly design with easy to see control panel, humidifier bottle holder and filter access.

Low noise level and neutral color help to blend into any home environment.

Reliable concentrator provides unlimited, worry-free refillable oxygen supply.

HCPCS Code: E1390*



1. Invacare® HomeFill® Oxygen System connection
2. HomeFill compressor and Ready Rack have been updated to coordinate with the Perfecto₂ Concentrator. Use with Invacare cylinders
3. Easy to maintain filter (not pictured)



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2019

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Lotus is a registered trademark of Inovo, Inc. Specifications are subject to change without notification.
Form No. 19-139
190908

| Model # | Dimensions | Weight | Shipping Weight | Quantity | Oxygen Flow | Flow Range | Maximum Outlet Pressure | Pressure Relief Mechanism |
|----------|--------------------------------|--|---------------------------|-----------------------------|---------------------------|-------------|-------------------------|---------------------------|
| IRC5PO2V | 23" H 15" W 11.5" D | 40 lb. | 45 lb. | 1 Per Case 27 Per Pallet | 5 LPM | 0.5-5 LPM | 5+/- 0.5 PSI | 35 PSI +/- 5 PSI |
| | | | | | | | | |
| | Oxygen Purity | Filter | Electrical Requirements | Rated Current Input | Average Power Consumption | Sound Level | Altitude | Limited Warranty |
| | 95.6%-87% At All Flow Rates | Cabinet Compressor Inlet Outlet HEPA | 120 VAC +/- 10%, 60 Hz | 3.0 A | 325 W @ 5 LPM | 45 db(A) | 8,000 ft. | 3 Years |

*All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.



Yes, you can.®