

ECLIPSE 5®

**TIME-PROVEN TRANSPORTABLE
OXYGEN DELIVERY.**



**Clinically effective.
Unmatched versatility.
Trusted performance.**


CAIRE®
O₂ is what we do.

SETTING THE STANDARD IN PORTABLE OXYGEN.

The award-winning Eclipse 5 transportable oxygen concentrator combines continuous and pulse dose flow options in one device, offering greater flexibility for hospital discharges, traveling and everyday portability. Its proven clinical performance allows the device to adapt to the patient, rather than the patient adapting to the device—making it the perfect companion for patients to stay active.

Clinically-Effective Oxygen Delivery



Smart technology to ensure saturation

autoSAT® Automatic adjustment to patient's respiratory rate

Adjustable Rise Time Maximizing inhaled oxygen

UltraSense® Sensitive breath detection technology

myCAIRE™ Telehealth solution



Continuous or pulse dose delivery to meet patients' needs

With 6 continuous flow settings up to 3 LPM and 9 pulse flow settings up to 192 mL bolus, the Eclipse 5 serves a wide range of patients.

Unmatched Versatility and Ease of Use



Compact and Transportable

Only 18.7 lbs, with an easy-to-attach traveling cart with telescoping handle and durable wheels.



Easy-to-Read Display

Large fonts and LCD display in LPM or mL for more accurate prescription verification.



Extended Battery Life

Over 5 hours of battery life with ability to recharge while plugged into DC power at most settings.



Enables Airline Travel

FAA-approved for commercial airline travel.



Quiet

At 40 dBA the sound level is one of the quietest on the market



24/7 Use

Meets a variety of needs for all phases of activity.

Proven Quality and Performance



Industry-Leading Reliability

Greater than 99% reliability rate out of the box.



Reduced Maintenance

Robust design features superior sieve bed technology reducing required maintenance.



Comprehensive Long-Term Warranty

CAIRE is the first manufacturer to offer a 3-year warranty on the device AND sieve beds.



Powerful Air Separation Technology

Vacuum Pressure Swing Adsorption (VPSA) is the most efficient way to deliver high purity oxygen.



PATIENT SATISFACTION BUILT-IN.

Large, easy-to-read, display



Customized oxygen delivery capabilities to meet broad range of prescriptions

Smart features to keep patients saturated

Long lasting battery life – enough to let patients stay active

Portfolio of accessories to get you anywhere you want to go

Removable cart for greater flexibility





Serving your patients across the continuum of CAIRE.

With innovative products that address the needs of new, highly ambulatory patients, later stage high-flow patients and everybody in between - all from one trusted source.



O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.

Learn more at:
eclipseoxygen.com



SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.


CAIRE[®]
O₂ is what we do.