



PHILIPS

Hospital
respiratory care

Confidence in your patient's comfort

Philips Respironics patient interfaces and circuits for noninvasive ventilation

The print quality of this copy is not an accurate representation of the original.



Philips Respironics AF531



AF531 with CapStrap



AF531 EE leak 1

Key benefits

- CapStrap headgear provides an excellent fit and simple reapplication
- Supports your infection control efforts with the CleanClip system
- You can adapt one mask to multiple ventilators with interchangeable elbows

Compatible with a leak 1 or leak 2 entrainment elbow (EE), EE NIVO elbow (EE/NE), standard elbow (SE), SE NIVO elbow (SE/NE), and bronchoscopy elbow (BE)
 For accessories, see page 23
 For headgear, see page 26

Philips Respironics AF531 with CapStrap single-use

Description	Quantity	Product number	Part number
AF531 single-use with EE leak 1 elbow, CapStrap headgear, small*	1	1061727	989805628611
AF531 single-use with EE leak 1 elbow, CapStrap headgear, medium	1	1061728	989805628621
AF531 single-use with EE leak 1 elbow, CapStrap headgear, large	1	1061729	989805628631
AF531 single-use with EE leak 1 elbow, CapStrap headgear, small*	10	1072627	989805628731
AF531 single-use with EE leak 1 elbow, CapStrap headgear, medium	10	1072628	989805628741
AF531 single-use with EE leak 1 elbow, CapStrap headgear, large	10	1072629	989805628751
AF531 single-use with SE standard elbow, CapStrap headgear, small*	1	1061730	989805628641
AF531 single-use with SE standard elbow, CapStrap headgear, medium	1	1061731	989805628651
AF531 single-use with SE standard elbow, CapStrap headgear, large	1	1061732	989805628661
AF531 single-use with SE standard elbow, CapStrap headgear, small*	10	1072630	989805628761
AF531 single-use with SE standard elbow, CapStrap headgear, medium	10	1072631	989805628771
AF531 single-use with SE standard elbow, CapStrap headgear, large	10	1072632	989805628781
AF531 single-use with EE leak 1/NE nebulizer elbow, CapStrap headgear, small*	10	1072618	989805635001
AF531 single-use with EE leak 1/NE nebulizer elbow, CapStrap headgear, medium	10	1072619	989805635011
AF531 single-use with EE leak 1/NE nebulizer elbow, CapStrap headgear, large	10	1072620	989805635021

* Intended for use by patients 7 years or older and > 40 lb (20 kg)



Philips Respironics AF531



AF531 with four-point headgear



AF531 Bronchoscopy*



AF531 SE

Key benefits

- Supports your infection control efforts with the CleanClip system
- You can adapt one mask to multiple ventilators with interchangeable elbows

For accessories, see page 23
For headgear, see page 26

Philips Respironics AF531 with four-point headgear single-use

Description	Quantity	Product number	Part number
AF531 single-use with EE leak 1 elbow, four-point headgear, small**	1	1061718	989805628521
AF531 single-use with EE leak 1 elbow, four-point headgear, medium	1	1061719	989805628531
AF531 single-use with EE leak 1 elbow, four-point headgear, large	1	1061720	989805628541
AF531 single-use with EE leak 1 elbow, four-point headgear, small**	10	1072621	989805628671
AF531 single-use with EE leak 1 elbow, four-point headgear, medium	10	1072622	989805628681
AF531 single-use with EE leak 1 elbow, four-point headgear, large	10	1072623	989805628691
AF531 single-use with SE standard elbow, four-point headgear, small**	1	1061721	989805628551
AF531 single-use with SE standard elbow, four-point headgear, medium	1	1061722	989805628561
AF531 single-use with SE standard elbow, four-point headgear, large	1	1061723	989805628571
AF531 single-use with SE standard elbow, four-point headgear, small**	10	1072624	989805628701
AF531 single-use with SE standard elbow, four-point headgear, medium	10	1072625	989805628711
AF531 single-use with SE standard elbow, four-point headgear, large	10	1072626	989805628721
AF531 single-use with EE leak 1 nebulizer elbow, four-point headgear, small**	10	1061633	989805634971
AF531 single-use with EE leak 1 nebulizer elbow, four-point headgear, medium	10	1061634	989805634981
AF531 single-use with EE leak 1 nebulizer elbow, four-point headgear, large	10	1061635	989805634991

* Bronchoscopy elbow sold separately, see page 23

** Intended for use by patients 7 years or older and > 40 lb (20 kg)



Philips Respironics AF531 with leak 2 elbow



Key benefits

- Leak built into the mask
- Can be used without an additional exhalation port
- Available with four point headgear and CapStrap headgear
- Convenient ten-packs

For accessories, see page 23
 For headgear, see page 26

AF531 EE leak 2

Philips Respironics AF531 with leak 2 elbow single-use

Description	Quantity	Product number	Part number
AF531 with EE leak 2 elbow, four-point headgear, small*	1	1061724	989805628581
AF531 with EE leak 2 elbow, four-point headgear, medium	1	1061725	989805628591
AF531 with EE leak 2 elbow, four-point headgear, large	1	1061726	989805628601
AF531 with EE leak 2 elbow, four-point headgear, small*	10	1100393	989805645651
AF531 with EE leak 2 elbow, four-point headgear, medium	10	1100394	989805645661
AF531 with EE leak 2 elbow, four-point headgear, large	10	1100395	989805645671
AF531 with EE leak 2 elbow, CapStrap headgear, small*	10	1100390	989805645621
AF531 with EE leak 2 elbow, CapStrap headgear, medium	10	1100391	989805645631
AF531 with EE leak 2 elbow, CapStrap headgear, large	10	1100392	989805645641

* Intended for use by patients 7 years or older and > 40 lb (20 kg)