

Breeze Nasal Comfort

The Breeze Comfort Range offers an innovative concept!
You can adjust it and then forget it!

	S	M	L
Masks*	M-268020-12	M-268030-12	M-268040-12
Cushion	M-429400-06	M-429400-07	M-429400-08
Forehead Support	M-329411-10		
Headgear	M-262040-10		
Headgear Clips	M-429440-01		

*The Nasal Comfort Mask is delivered in its cardboard packaging with 1 headgear fitted with its clips
1 storage bag
1 forehead support tab size
1 silicone cushion size.

Noise level : 27,5 dBA
Measured with an uncertainty of 2,5 dBA in accordance with standard EN ISO 17510-2:2009.



Breeze Pillows Mask

Thanks to minimal contact with the wearer's face, the Breeze Pillows Mask provides a real feeling of freedom whilst at the same time offering advantages that facilitate sleep, and greater freedom of movement .

	S	M	L
Masks*	M-268090-30		
Set of Prongs	M-329400-02	M-329400-03	M-329400-04
Headgear	M-262040-30		

* The Breeze Pillows Mask is delivered in its cardboard packaging with 1 headgear
1 storage bag
3 set of prongs sizes

Noise level : 22,5 dBA
Measured with an uncertainty of 2,5 dBA in accordance with standard EN ISO 17510-2:2009.

