



Taking comfort further

Philips Respironics AF811 oro-nasal mask with gel cushion

The Respironics AF811 oro-nasal single-use mask provides long-term comfort by combining a soft gel cushion with a pliable silicone seal. The gel cushion evenly distributes pressure on the face, and the silicone flap seals around the sensitive nasal bridge with minimal pressure. The CapStrap headgear provides additional stability and comfort while facilitating oral access.

Key advantages

- AF811 cushion uses a gel pad with a silicone flap for maximum comfort
- CapStrap headgear system simplifies mask application and routine patient care
- AF811 is available in three sizes to fit a wide patient population

PHILIPS

Philips Respironics AF811

The Philips Respironics AF811 oro-nasal single-use mask with gel cushion has additional features to improve patient comfort and noninvasive ventilation (NIV) performance for longer term patients in the institutional environment.

- Foam forehead pads provide patient comfort and mask stability.
- Swivel clips simplify re-application by eliminating the need to refit after removal.
- Adjustable forehead support improves overall mask fit and reduces skin breakdown.
- With a choice of the standard elbow (SE) for dual-limb systems and the entrainment elbow (EE) for single-limb systems, clinicians can use the AF811 mask with most noninvasive ventilators.
- Velcro® grab tabs allow quick, easy removal and adjustment of headgear.
- Integrated sizing gauge provides a convenient way to reduce cross-contamination and waste.
- Lower cost four-point replacement headgear extends mask life.



Entrainment elbow (EE)



Standard elbow (SE)

Description	Quantity	Part number
AF811 oro-nasal single-use mask EE, small	1	1063047
AF811 oro-nasal single-use mask EE, medium	1	1063016
AF811 oro-nasal single-use mask EE, large	1	1063015
AF811 oro-nasal single-use mask EE, small	5	1063053
AF811 oro-nasal single-use mask EE, medium	5	1063052
AF811 oro-nasal single-use mask EE, large	5	1063051
AF811 oro-nasal single-use mask SE, small	1	1063050
AF811 oro-nasal single-use mask SE, medium	1	1063049
AF811 oro-nasal single-use mask SE, large	1	1063048
AF811 oro-nasal single-use mask SE, small	5	1063056
AF811 oro-nasal single-use mask SE, medium	5	1063055
AF811 oro-nasal single-use mask SE, large	5	1063054
CapStrap replacement headgear	5	1019543
Four-point replacement headgear, petite	5	1018472
Four-point replacement headgear, small-medium-large	5	1007365
Nasogastric tube sealing pad	10	452040

Please visit www.philips.com/AF811



© 2009 Koninklijke Philips Electronics N.V.
All rights are reserved. Respironics and CapStrap are trademarks of Respironics, Inc. and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com
fax: +31 40 27 64 887

Printed in the U.S.A.
4522 962 54621 * SEP 2009

Philips Healthcare
Global Information Center
P.O. Box 1286
5602 BG Eindhoven
The Netherlands