

Invacare® Transport chair is designed for safe and efficient patient transport. The Transport chair has the added advantage of sharing interchangeable parts with the other chairs in the transport series. The Transport chair is compact and lightweight.

TRAN17FR and TRAN19FR Features

- Stylish, Silver Vein, powder-coated carbon steel frame
- Footrests are included with the chair and match frame finish
- Black nylon upholstery is flame retardant and bacteria resistant
- Full-length arms are padded for comfort
- Carry straps on seat allow for ease of folding and lifting
- Fold-down back makes transport and storage easy

TRAN19FR Transport Chair.



Specifications

Model numbers	TRAN17FR	TRAN19FR
Chair dimensions		
Seat width	17"	19"
Seat depth	16"	16"
Back height	19" fold-down	19" fold-down
Seat-to-floor height	19"	19"
Casters		
	8"	8"
Front riggings		
Swingaway footrest	17"-21" knee-to-heel	17"-21" knee-to-heel
Swingaway elevating legrest	17.5"-21" knee-to-heel	17.5"-21" knee-to-heel
Arm styles		
	Permanent full-length, padded	Permanent full-length, padded
Overall dimensions		
Overall width	19"	21"
Folded width	7.25"	7.75"
Overall depth w/o riggings	27"	27"
Overall height	39"	39"
Weight w/o riggings	24.5 lb.	25 lb.
Weight capacity	250 lb.	250 lb.
Limited warranty	12 months	12 months

Note: All specifications are approximate